

Now easily perform second surgeries!

SIS Bone Profiler

Useful Simple Effective



Features

- Excellent bone removal from the fixture neck
- Multiple diameters available
- Compatible with all implant systems

Access implant platform

- Easily remove bone blocking access to the fixture connection


Safe and easy to use

- No need for a guide screw or fixture attachment
- Remove bone and create the appropriate outline for healing abutments, impression copings, and temporary and final prosthetics
- Remove bone from the implant shoulder without damage to the fixture

※ Use only if the bone interferes with the abutment's emergence profile
If the abutment height is above 4mm, refrain from bone profiling

Bone Profiler Line-up

Internal Type (Luna S)



External Type (Sola S)

